

# Testa CHALLENGE

2020

Technical Description



# Technical description


A start to finish Bioprocess – Producing and purifying a domain antibody fragment, dAb.

## Process outline:

- Production of dAb in *E.coli* in bioreactor
- Release of dAb to periplasm under heat treatment
- Clarification and concentration of dAb using cross flow filtration
- Purification of dAb using affinity chromatography

More info on antibody fragments,  
suggested reading:

<https://www.mdpi.com/2073-4468/4/3/259/htm>



# Process overview

Cultivation in 50 L  
bioreactor

Clarification of harvest  
– Cross flow filtration

Purification of product  
– Affinity  
Chromatography

# Timeline

Preparation of culture medium  
Prep and installation of bioreactor.  
Thawing of cell bank and pre-culture start

Day 1

Cultivation in XDR-50  
Preparation of solutions for filtration  
Installation and preparation of harvest filter

Day 2


Heat treatment in Bioreactor  
Transfer of cell broth to filtration line  
Clarification of harvest in ÄKTA Ready

Day 3

Preparations for purification, buffers and method scripting  
Purification of product using ÄKTA pilot and protein L resin in AxiChrom column

Day 4

# Upstream process flow chart



# Clarification process flow chart

## Harvest using Cross Flow Filtration/Hollow fibers

### Material to process:

- 8.5 L *E.coli* cultivation
- Conc.: ~2 g dAb/L

### Cross flow filtration:

- Hollow fiber cartridge
- 750 kD cut off
- 8.5 L *E.coli* cultivation is processed in 1 run with 0.28 m<sup>2</sup> filter area
- 10 L recirculation bag

### Clarification conditions:

#### Buffers:

- 20 mM citrate, 0.8 M NaCl, pH 5.0
- 0.5 M NaOH, 300 ppm NaOCl
- 0.1 M NaOH
- 2.5 x concentration
- 3 x wash
- Feed flow 8L/min
- Permeate flow 70 mL/min
- 15 L clarified dAb solution

ÄKTA Ready Flux, with flowkit



UFP-750-E-6A

### Clarification method:

Step	Solution	Time (min)	Volume (L)
Rinse+normalized water flux test	Water	35	16.0
Buffer flush	20 mM citrate, 0.8 M NaCl, pH 5.0	15	3.0
Concentration	<i>E.coli</i> cultivation	60	8.5
Diafiltration	20 mM citrate, 0.8 M NaCl, pH 5.0	180	10.2
Empty waste	<i>E. coli</i> cells	10	3.4
Waste flush	20 mM citrate, 0.8 M NaCl, pH 5.0	10	3.0
Flush	Water	10	9.0
CIP	0.5 M NaOH + 300 ppm NaOCl	75	9.0
Rinse+normalized water flux test	Water	35	16.0
Storage	0.1 M NaOH	10	3.0

Total run time: 7 hours

### Permeate:

- Sterile filter to a 20 L bag
- 15 L dAb solution for chromatography

### Hollow fiber cartridge UFP-750-E-6A

- 0.28 m<sup>2</sup> filter area
- 63.5 cm long
- 170 fibers
- 1 mm lumen inner diameter

# Downstream process flow chart

## Capture step using Affinity Chromatography/ Capto L resin

### Material to process:

- 15 L permeate from clarification
- Conc.: 1 g/L
- 15 g dAb to purify

### Chromatography:

- Capto L resin
- Binding capacity: 14 g/L resin
- 15 g dAb is processed in 3 runs with a 360 mL column

### Chromatography conditions:

- Buffers:
  - 20 mM citrate, 0.8 M NaCl, pH 5.0
  - 20 mM citrate, pH 5.0
  - 20 mM citrate, pH 2.8
  - 15 mM NaOH
- Residence time during sample application: 4 minutes

ÄKTA pilot 600



AxiChrom 50/300



### Chromatographic method:

Step	Solution	CV	Comment
Eq.	20 mM citrate, 0.8 M NaCl, pH 5.0	1	300 cm/h
Sample appl.	dAb sample	5 L	4 min res.time
Wash 1	20 mM citrate, 0.8 M NaCl, pH 5.0	5	
Wash 2	20 mM citrate, pH 5.0	1	
Elution	20 mM citrate, pH 2.8	5	Peak fractionation
Wash 3	20 mM citrate, 0.8 M NaCl, pH 5.0	1	
CIP	15 mM NaOH	2	250 cm/h
Re-eq.	20 mM citrate, 0.8 M NaCl, pH 5.0	5	300 cm/h

Total run time: ~120 min.

### Eluate:

- Neutralize eluate ...?
- Collection bottle/bag ?

### AxiChrom 50/300:

- 18.3 cm bed height
- 360 mL Capto L chromatography resin



# Potential Analytics (TBD)

Offline cell count and viability (Vi-Cell XR)

Offline pH and pCO<sub>2</sub> (ABL9)

Glycerol and Acetate concentration (Cedex Bio)

dAb concentration (ÄKTA Pure 25)



# Instruments and consumables

Testa Challenge

# XDR MO single-use stirred-tank fermenter



<https://www.cytivalifesciences.com/en/us/shop/cell-culture-and-fermentation/microbial-fermentors/microbial-fermentation-systems/xcellerex-xdr-mo-single-use-stirred-tank-fermentors-p-06352>

# ÄKTA Ready Flux, Ultra filtration HF cartridge

150 cm



UFP-750-E-6A



<https://www.cytivalifesciences.com/en/us/shop/bioprocessing-filtration/tangential-flow-filtration/filtration-systems/akta-readyflux-p-07377>

<https://www.cytivalifesciences.com/en/us/shop/bioprocessing-filtration/tangential-flow-filtration/hollow-fiber-cartridges/laboratory-and-pilot-scale-ultrafiltration-hollow-fiber-cartridge-p-00452>

# ÄKTA pilot, AxiChrom 50/300

ÄKTA pilot 600



<https://www.cytivalifesciences.com/en/us/solutions/bioprocessing/products-and-solutions/downstream-bioprocessing/akta-pilot-600>

AxiChrom 50/300



<https://www.cytivalifesciences.com/en/us/shop/chromatography/columns/process-columns/axichrom-50-to-200-mm-chromatography-columns-p-06215>

# Instruments

Instrument	Manufacturer and model	Comment
XDR-50 MO system	Cytiva	With baffles and filter condenser
XDM/XDUO 100 Mixer	Cytiva	With Heating/cooling jacket, load cells
Temperature control unit (TCU)	Polyscience	Heating/cooling capacity: 9 kW/ 1.5 HP
Peristaltic pump	Watson Marlow, 520	For feed & filtration of medium to the vessel
Pressure Monitor with disposable pressure sensors	SciLog	To monitor pressure during medium filtration <b>Optional</b>
Shaking Incubator	Infors, Ecotron	For cultivation of seed culture.
Moisture analyser	Mettler-Toledo, HB 43	To determine dry content of the culture
Spectrophotometer	–	For OD measurements
Centrifuge	Sorvall	Sample prep for dry content measurements. <b>Optional</b>
Nutrient and Metabolite analyzer	Cedex bio (Roche diagnostics)	Measurements of Glycerol and acetate in culture. <b>Optional</b>
pH-meter	–	For bioreactor probe calibration and solution preparation
Sterile tube sealer	Cytiva	To seal thermoplastic tubing (Sizes: 1/4 to 1 inch OD)
Sterile tube fuser	Cytiva	To fuse thermoplastic tubing (Sizes: 1/4 to 1 inch OD)
ÄKTA Ready Flux	Cytiva	Used for clarification of harvest
ÄKTA Pure 150	Cytiva	Affinity chromatography using Capto L resin
ÄKTA Pure 25 or HPLC	Cytiva	For analysis of dAB titer. <b>Optional</b>
HiScale 16/40 column	Cytiva	Affinity chromatography using Capto L resin

# Consumables

Consumable part	Cat. No.	Amount	Comment
ULTA Pure HC filter (5 inch single-use NFF, 0.2 µm)	12410095	2 pc	For sterile filtration of culture medium and clarified dAb.
Disposable 10 L hanging bag for recirculation	12410222	1	Bag needed for the ÄKTA readyflux system
Disposable 20 L hanging bag for permeate	12410224	2	Sterile bag for collecting product
ReadyMate DACs hose barb	28936687	1 pack	Disposable Aseptic Connectors 3/8" Barb, for aseptic connections.
ReadyMate DACs with Mini-TC	28936695	1 pack	Disposable aseptic connector converting ReadyMate to TC.
pH probe standard	817-00144	1 pc	For XDR-50 bioreactor; can be used up to 12 autoclavation cycles.
DO probe standard	817-00145	1 pc	For XDR-50 bioreactor; can be used up to 12 autoclavation cycles.
Probe sheaths	888-0138	1 pack	Four probe sheaths pcs included in one package.
XDM-100 BASIC	888-0164-C	1 pc	Mixer bag for medium preparation in XDUO-100.
XDR-50 Microbial BBA	888-0235	1 pc	50 L Single-Use microbial bioreactor bag
Capto L, 25 mL	17547801	2 pcs	Affinity resin for Dab capture
Pilot Scale Ultrafiltration Cartridge	UFP-750-E-6A	1 pc	Hollow fiber cartridge for clarification of harvest
Flow Kit plus RM, ÄKTA readyflux	29187382	1 pc	Flow kit for ÄKTA ready Flux